

## EXPERIMENT 9

### THE SOLUBILITY OF CALCIUM SULFATE


In this experiment you will determine the solubility of calcium sulfate in water and in a second solution that contains 0.050 F  $\text{KNO}_3$ , both at room temperature. We will call these the A solution and the B solution. I made these as follows:

the A solution: Excess  $\text{CaSO}_4 \cdot 2\text{H}_2\text{O}$  was added to a jug of distilled water and mixed and shaken for a period of months. Solid  $\text{CaSO}_4 \cdot 2\text{H}_2\text{O}$  remains in the mixture. We may be reasonably sure that the solution is saturated in calcium sulfate.

the B solution: Excess  $\text{CaSO}_4 \cdot 2\text{H}_2\text{O}$  was added to a jug of 0.050 F  $\text{KNO}_3$  (in water) and mixed and shaken as above. Again, solid  $\text{CaSO}_4 \cdot 2\text{H}_2\text{O}$  remains. This solution is saturated in calcium sulfate. Note that the solution is not saturated in  $\text{KNO}_3$  but simply contains 0.050 F  $\text{KNO}_3$ . (If some solid  $\text{KNO}_3$  were added it would dissolve.)

The "solubility" of a substance is the concentration of a saturated solution of the substance. We have two saturated calcium sulfate solutions. Determination of the solubilities amounts to determination of the concentration of calcium sulfate in the saturated solutions. Because calcium sulfate supplies one mol of  $\text{Ca}^{2+}$  ions to a solution for each mole of calcium sulfate that dissolves, the solubilities of calcium sulfate are just the concentration of  $\text{Ca}^{2+}$  ions in the saturated solutions.

The method you will use to determine the  $\text{Ca}^{2+}$  concentration in the two saturated  $\text{CaSO}_4$  solutions involves complexometric titration with EDTA (ethylenediaminetetraacetic acid). The reaction is symbolized by:



SOME NOTES: EDTA acts as a complexing agent for essentially all of the +2, +3, and +4 metal ions. It does not complex the alkali ( $\text{Li}^+$ ,  $\text{Na}^+$ , etc.) metal ions. When complexed by EDTA, the metal ion is engulfed as if surrounded by some many clawed beast. This sort of complex is termed a chelate and EDTA titrations are sometimes termed chelometric titrations.

EDTA is a tetraprotic acid; it may lose up to four  $\text{H}^+$  ions in aqueous solution. The term "EDTA" is not specific about the number of protons present in the substance. Instead this term simply refers to the substance without reference to the extent of protonation. To describe the extent of protonation we employ the symbol "Y" as follows.  $\text{Y}^{4-}$  is the 4- charged anion of EDTA,  $\text{HY}^{3-}$  is the monoprotonated (3- charged) anion of EDTA, etc. Most metal ion complexes

with EDTA involve  $Y^{4-}$  species. As examples, the  $Mg^{2+}$  and  $Ca^{2+}$  complexes with EDTA have the formulas  $MgY^{2-}$  and  $CaY^{2-}$ , respectively.

The titration method is a simple one involving titration of a known volume of  $Ca^{2+}$  solution with EDTA reagent of known concentration. In order for this method to work (that is, to give a sharp reproducible endpoint very close to the equivalence point) it is necessary that (1) the reaction between  $Ca^{2+}$  and EDTA is essentially complete and (2) that we have a suitable indicator for the titration. We fulfill the first requirement by appropriately adjusting the pH of the titration mixture. Because  $Y^{4-}$  is the salt of a weak tetraprotic acid  $H_4Y$ , the extent of its complexation with metal ions depends on the pH. By adjusting the pH to about 10 we assure that the reaction of EDTA with  $Ca^{2+}$  is essentially complete. (Under more acidic conditions  $H^+$  ions would successfully compete with  $Ca^{2+}$  for  $Y^{4-}$ . At a pH of say 2-3 only a very small fraction of  $Ca^{2+}$  would be complexed by EDTA.) The second condition, that of a suitable indicator, is more difficult. Some indicators for  $Ca^{2+}$  are available and suitable for use at pH 10 but they tend to react rather slowly and are hard to use. To make matters worse, the color changes are hard to see. We will employ a trick called a displacement titration.

Consider a solution containing  $Ca^{2+}$  as the only species capable of complexing EDTA. Into the solution is added a certain amount of  $MgY^{2-}$  complex. The  $CaY^{2-}$  complex is stronger (has a larger formation constant) than the  $MgY^{2-}$  complex. As a consequence most of the  $Mg^{2+}$  is displaced from the complex. In other words, the reaction:  $Ca^{2+} + MgY^{2-} = CaY^{2-} + Mg^{2+}$  is a favorable one.

NOTE: The equilibrium constant for this reaction is readily calculated from known formation constants for  $Ca^{2+} + Y^{4-} = CaY^{2-}$  and  $Mg^{2+} + Y^{4-} = MgY^{2-}$  which values are  $4.9 \times 10^{10}$  and  $6.2 \times 10^8$ , respectively. The equilibrium constant for the displacement reaction is the ratio of these values:  $0.79 \times 10^2$ .

The analyte solution now contains  $Mg^{2+}$  ions which will react with added EDTA reagent. The endpoint of the  $Mg^{2+}$ , EDTA reaction is readily detectable by any one of a number of so-called "metallochromic" indicators. We will use one called Erio T (also called Eriochrome black T). The uncomplexed dye has a blue color at pH 10 but Erio T forms a red complex with  $Mg^{2+}$  (and many other 2+ and 3+ ions).

The titration proceeds as follows. First, the  $Ca^{2+}$  solution is adjusted to pH 10 and treated with  $MgY^{2-}$  solution. The displacement reaction releases most of the  $Mg^{2+}$  into solution. Erio T is added and the red color of the  $Mg^{2+}$ , Erio T complex is visible.

Now the titration is begun as EDTA reagent is added. The effect is to complex the now

free  $\text{Mg}^{2+}$  ions. Eventually enough EDTA is added so that both  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$  are present almost entirely as their complexes  $\text{CaY}^{2-}$  and  $\text{MgY}^{2-}$ .  $\text{Mg}^{2+}$  is no longer available for complexation with Erio T so the red color disappears and the blue color of free Erio T appears. Thus, the actual titration reaction is one between EDTA and  $\text{Mg}^{2+}$  with release of Erio T signaling the endpoint. However, the quantity of EDTA required to reach the endpoint depends only on the amount of  $\text{Ca}^{2+}$  initially present. In other words, we begin the titration with  $\text{MgY}^{2-}$  and  $\text{Ca}^{2+}$ . At the endpoint the solution contains  $\text{MgY}^{2-}$  and  $\text{CaY}^{2-}$ . The number of moles of EDTA required to cause this transformation is equal to the number of moles of  $\text{Ca}^{2+}$  present at the start.

We note at this point that any other substance present in the analyte solution that displaces  $\text{Mg}^{2+}$  from the  $\text{MgY}^{2-}$  complex will be titrated along with  $\text{Ca}^{2+}$ . As an example, suppose that a solution contains  $\text{Zn}^{2+}$  along with  $\text{Ca}^{2+}$ . The EDTA complex with  $\text{Zn}^{2+}$  is much stronger than the  $\text{MgY}^{2-}$  complex so that both  $\text{Ca}^{2+}$  and  $\text{Zn}^{2+}$  will displace  $\text{Mg}^{2+}$  from its complex. The quantity of EDTA reagent required to reach the endpoint will be the sum of the number of moles of  $\text{Zn}^{2+}$  and  $\text{Ca}^{2+}$  present. In our experiment this does not present a problem because  $\text{Ca}^{2+}$  is the only species available for reaction in the saturated calcium sulfate solutions.

### **BEFORE THE LABORATORY**

Read through the lab directions carefully. Make up a DATA SHEET with a very brief description of just what things need to be measured and how. Record your data on the SHEET and turn it in with your lab report.

## IN THE LABORATORY

### 1. Preparation of $\text{MgY}^{2-}$ solution.

It is not necessary to know the exact concentration of this solution but it is essential that  $\text{Mg}^{2+}$  and EDTA be present in a one - to - one mole ratio.

In two small labeled beakers take about 40-50 mL each of approximate 0.05  $\underline{\text{F}}$   $\text{MgSO}_4$  and 0.05  $\underline{\text{F}}$   $\text{Na}_2\text{H}_2\text{Y}$  from the laboratory stocks. (Beware. Two different concentrations of  $\text{Na}_2\text{H}_2\text{Y}$  solution are available.)

Measure about 30 mL of each solution into a small clean erlenmeyer flask. Add 1 drop of phenolphthalein to the flask and dropwise add 1  $\underline{\text{F}}$   $\text{NaOH}$  to the first pink color. (The  $\text{NaOH}$  is required for the complexation reaction,  $\text{H}_2\text{Y}^{2-} + \text{Mg}^{2+} + 2\text{OH}^- = \text{MgY}^{2-} + 2\text{H}_2\text{O}$ .) About 2-3 mL of  $\text{NaOH}$  solution will be needed. Now add about 10 mL of pH 10 buffer and 2-3 drops of Erio T indicator solution. At this point the solution may be red, which signals the presence of uncomplexed  $\text{Mg}^{2+}$ . The solution may be blue, which signals an excess of EDTA. If a miracle happens, the solution will be purple, signaling an essentially exact one - to - one mole ratio of  $\text{Mg}^{2+}$  to EDTA. (Red + blue = purple.) In any case your objective is to make a solution that will turn blue upon addition of 1 - 2 drops of EDTA solution and will turn red upon addition of 1 - 2 drops of  $\text{Mg}^{2+}$  solution. Depending on the color of your solution add dropwise and carefully and with adequate mixing either  $\text{Mg}^{2+}$  or EDTA until the color just turns. (THINK FIRST.) Test the solution by adding 1-2 drops of the other reactant to see that the color changes back. At this point you have the required  $\text{MgY}^{2-}$  solution. Label and cover the flask.

### 2. Standardization of 0.01 $\underline{\text{F}}$ $\text{Na}_2\text{H}_2\text{Y}$

Take about 100 mL of standard 0.01  $\underline{\text{F}}$   $\text{MgSO}_4$  solution (Record the actual concentration.) and pipet 20 mL portions into each of 3 erlenmeyer flasks. Add about 20 mL of water, 5 mL of  $\text{MgY}^{2-}$  solution, 3 mL of pH 10 buffer and 2-3 drops of Erio T to each flask. Rinse and fill a 25 mL buret with 0.01  $\underline{\text{F}}$   $\text{Na}_2\text{H}_2\text{Y}$  solution (NOT 0.05  $\underline{\text{F}}$   $\text{Na}_2\text{H}_2\text{Y}$  solution). Titrate to the blue color of Erio T. You should have three closely matched endpoint volumes,  $V_1^*$ . If not, make a few more titrations.

NOTE: Because you are using a  $\text{Mg}^{2+}$  solution to standardize the EDTA it might not seem necessary to add the  $\text{MgY}^{2-}$  solution. From the point of view of getting a good color change, the  $\text{MgY}^{2-}$  is not necessary. However, the conditions met in the standardization and the sample titrations should always be matched as closely as possible. Therefore, it is most wise to add all of the

same components here that will later be present in the sample titrations, including the  $\text{MgY}^{2-}$  solution.

### 3. Analysis of the Saturated $\text{CaSO}_4$ Solutions

At this point you have standardized the EDTA reagent and will use it to analyze the two saturated  $\text{CaSO}_4$  mixtures, A and B. These analyses will be made by the same titration procedure as above. That is, mixtures of sample ( $\text{CaSO}_4$  solution), water, buffer and  $\text{MgY}^{2-}$  solution are titrated with the EDTA reagent. Before doing so it is necessary to decide on how much of the saturated  $\text{CaSO}_4$  solutions should be measured out.

The solubility of  $\text{CaSO}_4$  (F.W. = 136.2 g/mol) in water is roughly 2 g/L near room temperature.

You will use a 25 mL buret for the titrations. In order to minimize the buret reading error, the sample size should be taken to give at least a 10 mL titration, but necessarily less than a 25 mL titration. In short you would like to use roughly 15 mL of EDTA reagent in a titration. Decide on an appropriate sample volume for the titrations. For the sake of simplicity use the same sample volume to analyze both the A and B solutions. Call this volume  $V_S$ .

a. Take an amount of saturated  $\text{CaSO}_4$  in water, A, from the laboratory stock appropriate to 3 or 4 titrations. **Be very careful to not shake or disturb the stock solution in any way!** Pipet  $V_S$  mL of the solution into titration flasks and add water, buffer, etc. as before. You should have 3 closely spaced endpoint volumes,  $V_2^*$ .

b. Take an appropriate amount of the saturated  $\text{CaSO}_4$  solution made up in 0.050  $\text{F KNO}_3$ , B, and measure  $V_S$  mL into three titration flasks. **Be very careful to not shake or disturb the stock solution in any way!** Titrate as before. You should have 3 closely spaced endpoint volumes,  $V_3^*$ .

## THE LABORATORY REPORT

1. Calculate the average  $V_1^*$ , its standard deviation,  $s_e$ , its standard error,  $s_m$ , and its 95 % CI for the standardization titrations.
2. Calculate the concentration of the EDTA titrant (0.01  $\underline{F}$   $\text{Na}_2\text{H}_2\text{Y}$ ). Determine the error in the [EDTA] by propagating the errors in the calculations based on the 95 % CI for  $V_1^*$  and the tolerances of the measuring devices (see Harris pgs 28-34 and see the opening preamble on the Guide to Statistics for how to propagate error).
3. Calculate the average  $V_2^*$ , its standard deviation,  $s_e$ , its standard error,  $s_m$ , and its 95 % CI for the titrations of the saturated  $\text{CaSO}_4$  in water.
4. Calculate the solubility of  $\text{CaSO}_4$  in water. Determine the error in the solubility of  $\text{CaSO}_4$  in water by propagating the errors in the calculations based on the 95 % CI for  $V_2^*$  and the tolerances of the measuring devices (see Harris pgs 28-34 and see the opening preamble on the Guide to Statistics for how to propagate error).
5. Calculate the average  $V_3^*$ , its standard deviation,  $s_e$ , its standard error,  $s_m$ , and its 95 % CI for the titrations of the saturated  $\text{CaSO}_4$  in 0.050  $\underline{F}$   $\text{KNO}_3$ .
6. Calculate the solubility of  $\text{CaSO}_4$  in 0.050  $\underline{F}$   $\text{KNO}_3$ . Determine the error in the solubility of  $\text{CaSO}_4$  in 0.050  $\underline{F}$   $\text{KNO}_3$  by propagating the errors in the calculations based on the 95 % CI for  $V_3^*$  and the tolerances of the measuring devices (see Harris pgs 28-34 and see the opening preamble on the Guide to Statistics for how to propagate error).
7. The object of this question is to calculate the "thermodynamic solubility product" for  $\text{CaSO}_4$ . We need to talk.

The thermodynamic solubility product for  $\text{CaSO}_4$  is symbolized  $K_{\text{SP}}$  (sometimes  $K_{\text{SP}}^\circ$ ) and defined by  $K_{\text{SP}} = (\mathbf{a}_{\text{Ca}})(\mathbf{a}_{\text{SO}_4})$ . The  $\mathbf{a}$ 's symbolize **activities** and are related to concentrations by  $\mathbf{a}_{\text{Ca}} = \gamma_{\text{Ca}}[\text{Ca}^{2+}]$  and  $\mathbf{a}_{\text{SO}_4} = \gamma_{\text{SO}_4}[\text{SO}_4^{2-}]$ . The  $\gamma$ 's are called the "activity coefficients" of the ions. In general the value of  $\gamma$  for any ion is not equal to 1.00 so that activities are not the same as concentrations. However, in very, very dilute solutions the values of ionic activity coefficients all approach 1.00. Activity coefficient values deviate from 1.00 at higher ionic

concentrations because the presence of other ions in the solution provides electrical attraction and repulsion interactions with any given ion. The effect of these interactions is to "shield" a given ion and so, in general, cause its apparent concentration, the activity, to differ from the actual concentration.

The **Debye-Huckel** theory allows calculation of activity coefficient values with good accuracy at ionic strengths below about 0.2 M. (And what, pray tell, is an ionic strength?)

The ionic strength of a solution is a numerical measure of the number and strength of electrical interactions that occur between ions in a given solution. Its value is symbolized by  $\mu$  or  $I$  and is defined by  $\mu = 1/2 \sum C_i z_i^2$ .  $C_i$  is the concentration of a given ion in a solution and  $z_i$  is the charge of that ion. The sum involves all of the ions in the solution.

AN EXAMPLE: A certain solution contains 0.0100 F  $\text{CaSO}_4$  in pure water. The only ions present in appreciable concentrations here are  $\text{Ca}^{2+}$  and  $\text{SO}_4^{2-}$  with concentrations of 0.0100 M each. The ionic strength is  $1/2[(0.0100)(+2)^2 + (0.0100)(-2)^2] = 0.0400$  M.

ANOTHER EXAMPLE: A certain solution contains 0.0100 F  $\text{CaSO}_4$  and 0.050 F  $\text{KNO}_3$  in water. The ions present are: 0.0100 M  $\text{Ca}^{2+}$ , 0.0100 M  $\text{SO}_4^{2-}$ , 0.050 M  $\text{K}^+$ , and 0.050 M  $\text{NO}_3^-$ . The ionic strength is  $1/2[(0.0100)(+2)^2 + (0.0100)(-2)^2 + (0.050)(+1)^2 + (0.050)(-1)^2] = 0.090$  M.

The activity coefficient for any given ion in the solution may now be calculated directly from the equation:

$$-\log \gamma_i = 0.511 z_i^2 \mu^{1/2} / (1 + 3.3 \times 10^{-3} \alpha \mu^{1/2})$$

$\gamma_i$  is the activity coefficient of a certain ion;  $z_i$  is the charge of that ion and  $\alpha$  is called the "ion-size" parameter and is essentially the size of the ion in units of nm. The numerical coefficients pertain to the solvent water and while they are slightly temperature-dependent the values given (0.511 and 3.3) are appropriate to the temperature range from about 20°C to about 30°C. Values of  $\alpha$  for  $\text{Ca}^{2+}$  and  $\text{SO}_4^{2-}$  are 0.6 nm and 0.35 nm, respectively.

Once the activity coefficient of an ion is found its **activity** is readily calculated from  $\mathbf{a}_i = \gamma_i[\text{conc}]$ .

AN EXAMPLE: What is the activity of  $\text{Ca}^{2+}$  ions in a solution of 0.0100 F  $\text{CaSO}_4$  in water?

We have already found  $\mu = 0.0400$  M. Substituting for  $\mu$  ( $= 0.0400$ ),  $z_i$  ( $= +2$ ) and  $\alpha$  ( $= 0.6$ ) into the Debye-Huckel equation gives  $-\log \gamma_{\text{Ca}} = 0.511(+2)^2(0.0400)^{1/2} / (1 + 3.3(0.6)(0.0400)^{1/2})$ ;  $-\log \gamma_{\text{Ca}} = 0.293$ ;  $\gamma_{\text{Ca}} = 0.510$

Now, the activity of calcium ions,  $\mathbf{a}_{\text{Ca}} = \gamma_{\text{Ca}}[\text{Ca}^{2+}]$

$$a_{\text{Ca}} = (0.510)(0.0100) = 0.00510.$$

The result in this example indicates that although the actual concentration of  $\text{Ca}^{2+}$  ions in the solution is 0.0100 M, it behaves - interacts with other things in the solution - as if it had a much lower concentration, 0.00510 M.

ANOTHER EXAMPLE: What is the activity of  $\text{Ca}^{2+}$  ions in a solution containing 0.0100 F  $\text{CaSO}_4$  in 0.050 F  $\text{KNO}_3$ ?

We already have found the ionic strength of this solution above;  $\mu = 0.090$  M. Substituting into the Debye-Huckel equation as before gives  $-\log\gamma_{\text{Ca}} = 0.385$ ;  $\gamma_{\text{Ca}} = 0.412$ . The activity of  $\text{Ca}^{2+}$ ,  $a_{\text{Ca}}$ , is  $(0.412)(0.0100) = 0.00412$ .

Consider the two examples. The solutions contain the same analytical concentrations of  $\text{Ca}^{2+}$  ions, 0.0100 M. However, the activity of  $\text{Ca}^{2+}$  ions in the second, higher ionic strength mixture, is lower than in the first. The presence of other ions in the solution influences the behavior of the  $\text{Ca}^{2+}$  ions.

Activity coefficients calculated from the Debye-Huckel equation are generally accurate to about  $\pm 1\%$  up to  $\mu = 0.1$  M for +1 and -1 ions. They are somewhat less accurate for divalent ions like  $\text{Ca}^{2+}$  and  $\text{SO}_4^{2-}$ , typically  $\pm 5\%$  at  $\mu = 0.1$  M. Furthermore, the  $\gamma$  values get worse and worse as the ionic strengths increase. At present there is no way to calculate reliable activity coefficient values in concentrated ionic solutions.

FINALLY, the question: Calculate the activity coefficients of  $\text{Ca}^{2+}$  and  $\text{SO}_4^{2-}$  in the A solution ( $\text{CaSO}_4$  in water) and use these to calculate the thermodynamic solubility product of  $\text{CaSO}_4$ . Your text gives a value for  $K_{\text{SP}}$ ; look it up. Did they get close?

8. Calculate the activity coefficients of  $\text{Ca}^{2+}$  and  $\text{SO}_4^{2-}$  for the B solution ( $\text{CaSO}_4$  in 0.050 F  $\text{KNO}_3$ ) and calculate a value for  $K_{\text{sp}}$ .

We wish to compare this value with the one in 4, and recognize that both values depend on calculated activity coefficients. We expect these ( $\gamma_{\text{Ca}}$  and  $\gamma_{\text{SO}_4}$ ) to be accurate to about  $\pm 5\%$  and will be satisfied if the two  $K_{\text{sp}}$  values match to within about  $\pm 10\%$ . Is this the case?

9. (EXTRA CREDIT - 3 points) A certain solution is prepared by adding together 10.0 mL of 0.0200 F  $\text{Ca}^{2+}$ , 5.00 mL of 0.0200 F  $\text{MgY}^{2-}$ , pH 10 buffer and water to make 50.0 mL.



(This solution is suspiciously like ones you have titrated.)  $\text{Ca}^{2+}$  displaces some  $\text{Mg}^{2+}$  from the EDTA complex. What percentage of  $\text{Mg}^{2+}$  remains as  $\text{MgY}^{2-}$ ?

HINT: Essentially all of the 0.100 mmol of EDTA present exists in either  $\text{MgY}^{2-}$  or  $\text{CaY}^{2-}$  complexes. The total concentration of EDTA is  $0.100 \text{ mmol}/50.0 \text{ mL} = 2.00 \times 10^{-3} \text{ F}$ . Consequently  $[\text{MgY}^{2-}] + [\text{CaY}^{2-}] = 2.00 \times 10^{-3}$ . Use this and mass balance equations for  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$  along with the equilibrium constant  $\text{Ca}^{2+} + \text{MgY}^{2-} = \text{CaY}^{2-} + \text{Mg}^{2+}$   $K = 0.79 \times 10^2$

to solve for  $[\text{MgY}^{2-}]$ .

10. (EXTRA CREDIT - 2 points) The pH is defined as  $-\log a_{\text{H}}$  where  $a_{\text{H}}$  is the activity of  $\text{H}^+$  ions.
- What is the pH value of an aqueous solution containing  $1.00 \times 10^{-2} \text{ F}$  HCl?
  - What is the pH value of an aqueous solution containing  $1.00 \times 10^{-2} \text{ F}$  HCl and  $0.10 \text{ F}$  KCl?

Note to instructor:

Preparation of the MgEDTA solution is problematic. It is important that there is no excess EDTA of  $\text{Mg}^{2+}$ . However, achieving this in a precise manner by adjusting to a purple color is not possible. The experiment could be improved in one of two ways.

- Buy a MgEDTA reagent instead of preparing your own.
- Instead of erichrome black T and the titration of  $\text{Mg}^{2+}$  at  $\text{pH} = 10$ , you could perform the titration of  $\text{Ca}^{2+}$  directly at  $\text{pH} 12$  using hydroxynaphthal blue (I think).

Also be sure to warn the students not to shake or disturb the stock solutions.