


Magnetic Moment

- In this experiment you will determine the magnetic susceptibility of one of the compounds you previously synthesized; e.g., $[\text{Cr}(\text{NH}_3)_6](\text{NO}_3)_3$, $\text{Mn}(\text{acac})_3$.
- From the corrected molar susceptibility you will calculate the magnetic moment of the compound, μ , which you will compare with the ideal spin-only moment for the compound.
- From $\mu_s = \sqrt{n(n+2)}$ B.M., the expected values of μ_s for transition-metal complexes with $n = 1-5$ unpaired electrons are as follows.

n	S	μ_s (B.M.)
1	1/2	1.73
2	1	2.83
3	3/2	3.87
4	2	4.90
5	5/2	5.92

Magnetic Susceptibility with an Evans Balance

- This experiment uses a modified form of the Guoy balance method, using a microscale apparatus devised by D. F. Evans and manufactured by Johnson-Matthey.



- A moveable magnet attached to a torsion balance detects the force created by diamagnetic and paramagnetic moments in the sample.
 - ▼ Diamagnetic moment makes the magnet move down.
 - ▲ Paramagnetic moment makes the magnet move up.

Evans Balance Susceptibility Calculation

- The gram susceptibility of the sample is given by

$$\chi_g = \frac{L}{m} \left[C(R - R_o) + \chi_v' A \right]$$

where χ_g = mass magnetic susceptibility

L = sample length in centimeters

m = sample mass in grams

C = balance calibration constant

R = reading from the digital display when the sample (in the sample tube) is in place in the balance

R_o = reading from the digital display when the empty sample tube is in place in the balance

χ_v' = volume susceptibility of air (0.029×10^{-6} erg·G⁻²·cm⁻³)

A = cross-sectional area of the sample

- The volume susceptibility of air can be ignored with solid samples, so the equation becomes

$$\chi_g = \frac{CL(R - R_o)}{m \times 10^9}$$

- ☛ In this experiment you will determine the value of C by taking readings on a sample of $[\text{Ni}(\text{en})_3]\text{S}_2\text{O}_3$, for which $\chi_g = 1.104 \times 10^{-5}$ erg·G⁻²·cm⁻³.

Calculating Magnetic Moment, μ

- The *molar magnetic susceptibility*, χ_M , is obtained from the mass magnetic susceptibility by multiplying by the molecular weight of the sample in units of g/mol; i.e.,

$$\chi_M = M\chi_g$$

- Units of χ_M are $\text{erg}\cdot\text{G}^{-2}$.
- Diamagnetic corrections need to be applied to this measured molar magnetic susceptibility.
 - The diamagnetic contributions arise from core paired electrons, ligand electron pairs, and counter ion electron pairs.

$$\chi_A = \chi_M - \{\chi_M(\text{core}) + \chi_M(\text{ligand}) + \chi_M(\text{ion})\}$$

- The diamagnetic correction factors are tabulated values, called Pascal's constants.
- ☛ Use the newly published values of G. A. Bain and J. F. Berry, *J. Chem. Educ.*, **2008**, 85, 532.
- We will assume that paramagnetic coupling is minimal ($\theta = 0$), so the simpler Curie Law equation applies:

$$\mu = 2.828\sqrt{\chi_A T}$$

- ☛ Do not forget to take the temperature at the time of the measurements to use in this calculation.

Making the Measurements Obtaining the Instrument Constant, C

- ① Determine R_o .
 - Zero the balance.
 - Weigh a clean, dry, empty sample tube on an analytical balance.
 - Place the tube in the Evans balance and read R_o .
 - The needle on the scale should drift no more than ± 1 .

- ② Fill the tube to at least 1.5 cm with $[\text{Ni}(\text{en})_3]\text{S}_2\text{O}_3$.
 - Gently tap the sample on a hard surface (away from the balance) to pack it well.
 - Measure the length of the sample in the tube in centimeters to obtain the value of L .
 - Obtain the mass of the tube and sample on the same analytical balance you previously used.
 - Rezero the balance, insert the sample, and take a reading to obtain R .
 - If a reading is off-scale, change to the X10 range and multiply the reading by 10.

- ③ Calculate the instrument constant, C , from your L , R_o , R , m data and the value of the mass susceptibility of $[\text{Ni}(\text{en})_3]\text{S}_2\text{O}_3$, $\chi_g = 1.104 \times 10^{-5} \text{ erg}\cdot\text{G}^{-2}\cdot\text{cm}^{-3}$, using

$$\chi_g = \frac{CL(R - R_o)}{m \times 10^9}$$

Making the Measurements Obtaining the Sample Data

- ① Empty the sample of $[\text{Ni}(\text{en})_3]\text{S}_2\text{O}_3$. (If absolutely sure it has not been contaminated, you may put it back in the original sample bottle.)
 - Clean and dry the sample tube.

- ② Redetermine R_0 .
 - Zero the balance.
 - Weigh the clean, dry, empty sample tube on an analytical balance.
 - Place the tube in the balance and read R_0 .
 - The needle on the scale should drift no more than ± 1 .

- ③ Fill the tube to at least 1.5 cm with your sample; e.g., $[\text{Cr}(\text{NH}_3)_6](\text{NO}_3)_3$, $\text{Mn}(\text{acac})_3$.
 - Gently tap the sample on a hard surface (away from the balance) to pack it well.
 - Measure the length of the sample in the tube in centimeters to obtain the value of L .
 - Obtain the mass of the tube and sample on the same analytical balance.
 - Rezero the balance, insert the sample, and take a reading to obtain R .
 - If a reading is off-scale, change to the X10 range and multiply the reading by 10.

- ④ Be sure to take the temperature near the balance to obtain T , needed for calculating μ by the Curie Law equation.

Calculating χ_g and μ for the Sample

- Calculate the mass susceptibility of your sample from L , R_o , R , m and the previously determined instrument constant, C , using

$$\chi_g = \frac{CL(R - R_o)}{m \times 10^9}$$

- Calculate the molar susceptibility by $\chi_M = M\chi_g$.
- Calculate the corrected molar susceptibility by subtracting the Pascal's constants:

$$\chi_A = \chi_M - \{ \chi_M(\text{core}) + \chi_M(\text{ligand}) + \chi_M(\text{ion}) \}$$

- Calculate μ from χ_A and the temperature in kelvin (K), using the Curie Law:

$$\mu = 2.828\sqrt{\chi_A T}$$